



PBI

Quick Reference Guide

QRS-Card™ Resting ECG System





Acquiring a New Resting ECG Test

1. Confirm that hardware and software are installed and working correctly (see separate installation instructions or user manual).
2. **Prep the patient.** *A good patient prep is most beneficial for reducing artifact and noise on the ECG
 - a. **Shave hair** if necessary from electrode sites.
 - b. **Rub skin with alcohol** for better electrode contact with skin.
 - c. **Let dry fully.**
3. **Attach electrodes to patient.** See electrode placement diagram below.
4. **Instruct the patient to relax.**



5. **Select the Resting ECG Icon**  **OR use hotkey F2.**
 - a. Depending on your setting (see user manual for setup instructions), the software will perform a number of pre-test tasks (may display setup screen or electrode position diagram).
6. Once you navigate through pre-test screens, **enter patient demographics in the 'Resting ECG Patient Information' card, click 'Ok'** to begin acquiring an ECG.
7. **A 10 Second ECG test will be acquired and stored in your ECG database.** If using default settings, for other settings see applicable User Manual.



Electrode Position Diagram


- V1 (Red) = Forth Intercostal space at the right border of the sternum
- V2 (Yellow) = Forth Intercostal space at the left border of the sternum
- V3 (Green) = Midway between locations V2 & V4
- V4 (Blue) = At the mid-clavicular line in the fifth intercostal space
- V5 (Orange) = At the anterior axillary line on the same horizontal level as V4
- V6 (Purple) = At the mid-axillary line on the same horizontal level as V4 & V5

Reviewing a Resting ECG Test

1. To view (or edit) a Resting ECG test, **highlight the desired test from the Cardiology Suite Master Screen**

test database and click  OR double click on the test to be reviewed.

2. The test record will be displayed in the Cardiology Suite Resting Dataview screen. The default view is a 4 x 3 lead view (see example below).



- Used to view or edit information about the patient.



- Shows various views of the data, including views of leads, complexes and rhythms. Also includes patient demographics, interpretation, ECG measurements and other test variables.



- provides summarized data for ECG measurements and automatic ECG Interpretation.

Also allows user to move measurement markers and edit interpretation statements.

Notice: QRS-Card Cardiology Suite can provide automatic computer interpretations. A computer generated interpretation cannot replace sound medical reasoning by a trained professional. Therefore, a physician should always review the interpretation and final analysis should always be obtained by a Physician. Proper administration, diagnosis and implementation of the test is the Physician's responsibility.



- allows user to print various Resting ECG reports to printer, PDF or Microsoft® Word.



QRS-Card™ Cardiology Suite Resting Quick Reference Guide

Resting, Stress , Holter ECG & Blood Pressure

Turn your Notebook, Tablet or Desktop PC into a Medical Device

Pulse Biomedical, Inc
112 Ivy Lane, King of Prussia, PA 19406
Phone: (610) 666-5510
Fax: (610) 666-5630

Pulse Biomedical GmbH
Wilhelm Bihler Str 4
72474 Winterlingen, Germany
Phone: (49) 7434-316038
Fax: (49) 7434-316039

www.QRSCard.com