



PBI


Quick Reference Guide

QRS-Card™ Stress ECG System



Acquiring a New Stress ECG Test

1. Confirm that hardware and software are installed and working correctly (see separate installation instructions or user manual).
2. **Prep the patient.** *A good patient prep is most beneficial for reducing artifact and noise on the ECG
 - a. **Shave hair** if necessary from electrode sites.
 - b. **Rub skin with alcohol** for better contact for electrodes.
 - c. **Let dry fully.**
3. **Attach electrodes to patient.** See electrode position diagram below.
4. **Instruct the patient to remain as relaxed as possible.**
5. **Select the Stress ECG Icon**  **OR use hotkey F3.**
 - a. Depending on your setting (see user manual for setup instructions), the software will perform a number of pre-test tasks (may display setup screen or electrode position diagram).
6. Once you navigate through these screens, **enter patient demographics in the 'Stress ECG Patient Information' card, click 'Ok'** to begin acquiring an ECG.
7. Stress ECG acquisition will follow protocol defined in Stress Setup screen. **Manually move through the pre-test phases** (Supine, Hyperventilation and Warm-Up) by **clicking the**  **advance stage button.**
8. **Exercise phase will begin** and continue according to protocol.
 - a. Other tasks such as **manually taking blood pressure, saving, annotating and printing ECG strips** can be executed during the ECG acquisition.
9. Once target heart rate is achieved or patient is fatigued, **click the**  **Begin Recovery button.**




Electrode Position Diagram

- V1 (Red)= Forth Intercostal space at the right border of the sternum
- V2 (Yellow)= Forth Intercostal space at the left border of the sternum
- V3 (Green) = Midway between locations V2 & V4
- V4 (Blue) = At the mid-clavicular line in the fifth intercostal space
- V5 (Orange) = At the anterior axillary line on the same horizontal level as V4
- V6 (Purple) = At the mid-axillary line on the same horizontal level as V4 & V5

Reviewing a Stress ECG Test

1. To view (or edit) a Stress ECG test, **highlight the desired test from the Cardiology Suite Master Screen**

test database and click  **OR double click on the test** to be reviewed.

2. The test record will be displayed in the Cardiology Suite Stress Dataview screen. The default view is a 4 x 3 lead view (see example below).



- Used to view or edit information about the patient.



- Shows various views of the data, including views of leads, complexes and selected strips.



- provides summarized data and tables for ECG measurements. Also allows user to edit interpretation statements.



- allows user to print various Stress ECG reports to printer, PDF or Microsoft® Word.



The Digital Cardiology Experts

Resting, Stress , Holter ECG & Blood Pressure
Turn your Notebook, Tablet or Desktop PC into a Medical Device

Pulse Biomedical, Inc
112 Ivy Lane, King of Prussia, PA 19406
Phone: (610) 666-5510
Fax: (610) 666-5630

Pulse Biomedical GmbH
Wilhelm Bihler Str 4
72474 Winterlingen, Germany
Phone: (49) 7434-316038
Fax: (49) 7434-316039

www.QRSCard.com